




DeTacMed

a brand of
DEHAS Medical Systems GmbH

Backboard Immobilization



MADE
IN
GERMANY

TACTICAL EMERGENCY MEDICINE | MILITARY | POLICE | RESCUE ORGANIZATIONS


WE THINK PRESSURE & FLOW

Spineboards

Our spineboards and basket stretcher are widely used for safely transporting seriously injured and wounded. They can float on water, be allowed to use during X-rays examination and also used together with head immobilizer. Made of high strength engineering plastic, the spine boards are enduring, anti-aging, weight loading and easy for storage. They are necessary tools in first-aid equipment for medical centers. Option to customize the color start from 1000 ea..


DTM-MS-BB1Y

Length: 183 cm
Width: 43 cm
Height: 5 cm
Material: HDPE
Color: yellow




DTM-MS-BB4

Length: 185 cm
Width: 46 cm
Height: 4 cm
Material: HDPE
Color: orange



DTM-MS-BB5Y


Length: 183 cm
Width: 46 cm
Height: 7,6 cm
Material: HDPE
Color: yellow



DTM-MS-BB7Y

Length: 185 cm
Width: 40 cm
Height: 5 cm
Material: HDPE
Color: yellow


*Pay load up to 300 kg,
14 hand-holds, 10 pins
weight 6,5 kg*



DTM-MS-BB9Y (CF)

Length: 183cm
Width: 40 cm
Height: 2,5 cm
Material: HDPE,
Carbon
fiber
Color: yellow


*Contoured design,
compact, rigid,
lightweight, extremely
durable, easy to carry
and store,
14 hand-holds, 12 pins,
weight 4.7 kg*



DTM-MS-PED

Length: 137 cm
Width: 35 cm
Height: 5 cm
Material: HDPE
Color: yellow

*Pediatric use
Oxford leather
stretcher surface,
light weight, small
size, easy to carry,
safe to use*



Basket Stretcher

DTM-MS-BS1000

Length: 218 cm
Width: 63 cm
Height: 18 cm
Material: ABS plastic,
stainless
steel
Color: orange


*Helicopter use,
special carrying strap
equipment, adjustable
foot safety, safety belts,
mattress, helicopter
fixation, fireproof
materials, no toxic or
environmentally harmful
substances, high wear
and corrosion resistance*



DTM-MS-BS2000

Length: 218 cm
Width: 63 cm
Height: 18 cm
Material: ABS plastic,
stainless
steel
Color: orange

*Split Type,
Helicopter use,
special carrying strap
equipment, adjustable
foot safety, safety belts,
mattress, helicopter
fixation, fireproof
materials, no toxic or
environmentally
harmful substances,
high wear and
corrosion resistance*





Scoop Stretcher

DTM-MS-SCP124

Length: 203 cm
Width: 42 cm
Height: 7 cm
Material: PE
Color: yellow



Plastic scoop stretcher
No rolling or lifting of the patient, immobilization in the position found, center can be opened, X-ray while secured on the stretcher, foldable for transport, head section recessed, ensure correct alignment of the cervical spine, the patient's jaw remains in the position found



DTM-MS-SCP223

Length: 168cm
Width: 44 cm
Height: 6,5 cm
Material: aluminium
Color: silver

Aluminium scoop stretcher
No rolling or lifting of the patient, immobilization in the position found, center can be opened, X-ray while secured on the stretcher, foldable for transport, separates in half during application and removal





Sheets

DTM-22011B

Length: 203 cm
Width: 101 cm
Material: enhanced fabric
vinyl, nylon
web supports
Color: black
Quantity: 10/case



14 handles transfer sheet,
easily foldable for compact storage in carry pouch, option to use carry handles or poles (not included)



DTM-22011W

Length: 203 cm
Width: 101 cm
Material: non-woven fabric 95 g
film 30 g
vinyl, nylon
web supports
Color: white
Quantity: 10/case

12 handles transfer sheet,
easily foldable for compact storage in carry pouch, option to use carry handles or poles (not included)



DTM-7033

Length: 228 cm
Width: 91 cm
Material: spundbond non-woven fabric
Color: blue
Quantity: 50/case

XPS fitted sheet
Designed specifically for the Stryker® XPS mattress, transfer of patients up to 180 kg



Head Immobilizer

Conforms to any type of backboard/scoop-stretcher, immobilizes the patient's neck in the position found, delivered complete with base, two support pads plus head-and chin-straps.

DTM-MS-IMMO-200A

Length:	38 cm
Width:	18 cm
Height:	25 cm
Material:	foam
Color:	red
Quantity:	8/case

Unique head strap designed, conforms patient's forehead, large ear holes for monitoring the patient, minimal interference with X-ray, MRI/CT, waterproof plastic coating, prevent bacterial growth, no blood or bodily fluids absorb, easy to clean and decontaminate, meets FMVSS 302 flammability rating




Strap Systems

DTM-MS-ST10

Length:	A 86 (± 5) cm
	B 147 (± 5) cm
	C 146 (± 5) cm
	D 126 (± 5) cm
	E 106 (± 5) cm
	F 123 (± 5) cm
Width:	G 5 (± 1) cm
Material:	PP fibre, Velcro
Color:	multi-colored
Quantity:	25/case

10 points of contact, excellent patient security and immobilization, high functionality, easy to use single piece straps



Rescue Knife

DTM-K-2686

Length:	x cm (open)
Material:	
Blade:	3cr13
Handle:	aluminium
Color:	black
Quantity:	120/case

Dedicated to police, paramedic and rescue professionals



20+ years of competence and solutions

Since 2002, we have specialized in the development and sale of medical technology for use inside and outside clinical facilities.

Information, contact or find dealers:



Oxygen Medical Products

Focused on medical oxygen, we develop and produce products for clinical and non-clinical use. Find out more:



DEHAS Medical Systems GmbH
Wesloer Straße 107-109
D-23568 Lübeck
Tel.: +49 451 80904-0
E-Mail: info@dehas.de
Web: www.DEHAS.de

WE THINK PRESSURE & FLOW